MESA

PRODUCER OF SPECIAL ALLOYS SINCE 1975

Rev.01-04/2018

Mesa di Sala Giacomo & C. S.n.C. via
dell'Artigianato, 35/37, 25039 Travagliato (BS)
ITALY Tel. +39 0306863251 Fax +39 0306863252
info@mesaitalia.it sales@mesaitalia.it www.mesaitalia.it

MAGNUM SOLARE CE 0123

LOT 00000

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MAGNUM SOLARE CE 0123

LOT 00000

Cobalt based dental alloy for ceramic, type 4*

*According to ISO 22674:2016

The alloy MAGNUM SOLARE is produced in conformity with standards ISO 9693:2012 and ISO 22674:2016, it is highly corrosion resistant, it does not contain Nickel and it is completely free of toxic elements like beryllium, indium and gallium.

Instruction for use

Modeling

Modeling with a minimum thickness of 0.3mm on simple structures and 0.5mm in complex cases or with patients with bruxism. Avoid elbows.

Manufacturing

Separate the milled structures from the disk with tungsten carbide burs or with separating discs, clean with ultrasonic and acetone or bioalcohol to eliminate any oil residues from the working process.

Continue manufacturing process using tungsten carbide burs.

The worked framework has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steam clean. After cleaning, the framework should not be touched anymore and should be held with hemostatic forceps only.

Oxidation

Fire for 5 min. under vacuum at 950-980°C. After firing, the oxide layer has to be carefully sandblasted using disposable oxide-aluminium equipment 110 to 150 µm at a pressure of 2.5 - 4 bar. Steam and check that the surface has an homogeneous grey surface.

Ceramization

It's possible to use Bonding. Apply the ceramic, following its manufacturer's instructions. Slow cooling is recommended.

Soldering

We suggest using our MAGNUM MAGNUM SALDATURA Co to weld.

Waste disposal

The processing scrap must be disposed of as special waste in accordance with the EC directives 2008/98/CEE on waste, and 94/62/CEE on packaging and packaging waste and in compliance with national legislation in force on the subject.

SAFETY INSTRUCTIONS

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Metal dusts and fumes are dangerous for health. Use an adequate exhaust fans during milling and polishing.
- Before prosthesis application verify if other metal implants are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- We identify every batch with a number. We recommend to write it down in patient's file to allow its complete traceability.
- This product does not need any special preservation precautions. MESA keeps batch traceability for 15 years.
- It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results.

Percentage chemical composition (m/m)	
Co	66
Cr	27
Others	Si, Mn

Physical and mechanical properties	
Density	8.4 g/cm ³
Solidus-liquidus temperature	1307-1417°C
Melting point	1470°C
Yield load strength (Rp 0.2)	395 MPa
Percentage elongation at fracture	11%
Modulus of elasticity	233 GPa
Vickers hardness	255 HV10
Colour	White
Thermal expansion coefficient 25-500°C	14.3 · 10 ⁻⁶ K ⁻¹
Thermal expansion coefficient 25-600°C	14.5 · 10 ⁻⁶ K ⁻¹
Maximum cooking temperature	980°C
Suggested ceramics	VITA VM13
Cytotoxicity test according to ISO 10993-5	Passed