RU ИЗДЕЛИЕ: STARBOND T14 DISC

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DIN EN ISO 14971: DIN EN ISO 22674: DIN EN ISO 5832-2: ASTM F 67

EN PRODUCT: STARBOND Ti4 DISC

Milling blank made of grade 4 titanium for CAM milling of fixed and removable and dental prostheses. According to the requirements as stipulated by ISO 22674 this alloy is a type 4 alloy. This alloy should be used by qualified and trained personnel for the designated scope of application.

Identified uses: removable, partial dentures, brackets, veneered crowns, bridges in the anterior region with up to three pontics, in the posterior region with up to two pontics; crosspieces, fastenings, implant-borne superstructures.

NOMINAL ANALYSIS IN MASS PERCENT:

Nominal values of alloy composition		
Ti	Other constituents:N; C; H; Fe; O	
>99 %	< 1 %	

ALLOY PROPERTIES:

Yield strength (Rp0.2)	504 MPa	Density	4.5 g/cm³
Tensile strength	599 MPa	Melting interval	1660°C
Elongation	23.5%	Thermal expansion coefficient 20-600°C	9,7 x 10 ⁻⁶ K ⁻¹
Vickers strength	>200 HV 5/30		

MODEL:

The crown walls of the virtual model should be contoured no thinner than 0.4 mm, so that the wall thickness after finishing and prior to ceramic or plastic veneering is at least 0.3 mm. Crowns and pontics are to be designed in accordance with the anatomical form of the teeth to provide for a consistent ceramic layer. Avoid formation of sharp ridges and undercut areas. Crosspieces between the pontics should be formed as strong and as high as possible. For the anterior region of bridges with up to three pontics, we recommend a minimum cross-section for the connectors of 5.5 mm². For the posterior region with up to two pontics, we recommend a minimum cross-section for the connectors of 8 mm².

GENERALLY

Frameworks can be finished using tungsten carbide burs suitable for use on titanium, according to the manufacturer's instructions for use.

CUTTING FRAMEWORKS

Frameworks can be cut using cross-cut tungsten carbide burs or cutters suitable for use on titanium.

FINISHING AND CLEANING

Finish frameworks with a clean tungsten carbide bur suitable for use on titanium. To prevent blistering caused by overlaps when firing the ceramic, tools should only be drawn unidirectionally over the surface. Also observe the manufacturer's recommended maximum speed for the instruments. Then blast the surface with pure aluminium oxide (110 - 250 µm) at a pressure of 2 to 3 bar. Finally, wash the framework thoroughly under running water or steam clean it. Degrease with ethyl alcohol, if required. Never use hydrofluoric acid!

VENEERING

Starbond Ti4 frameworks can be veneered with a veneering ceramic suitable for use on titanium. Please observe the corresponding ceramic manufacturer's instructions.

SIDE EFFECTS:

May be experienced in rare cases, such as allergies to alloy constituents or electrochemically-induced paresthesia.

INTERACTIONS:

None known.

CONTRAINDICATIONS:

For known intolerances, allergies to alloy constituents.

WARNING:

Biocompatible materials (ceramic) - please observe the ceramic manufacturer's instructions.

WRARRANTY:

Our recommendations for use whether given verbally, in writing or by practical instructions, are based on our own trials and experience and can only be considered as standard values. Our products are subject to constant further development. Therefore alterations in construction and composition are reserved.

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